




## Posterior-Inferior Retractors



The posterior-inferior retractor is placed with the point at 6 o'clock and the retractor's axilla resting on the ischium. The remaining blade of this retractor is used to retract the remaining capsule from the posterior lip of the acetabulum. Features a toothed inner arc for a better grip.



**PRODUCT NO'S:**

**7625-01** [Small Right]  
Overall Length: 10.75" (27.3 cm)  
Handle-to-Bend Length: 5.5" (14 cm)

**7625-02** [Small Left]  
Overall Length: 10.75" (27.3 cm)  
Handle-to-Bend Length: 5.5" (14 cm)



Toothed for a better grip



## Cobra Retractor with Hand Rest

The pointed cobra retractor with hand rest is placed anteriorly to help retract and expose.

**PRODUCT NO:**

**7130**  
Overall Length: 10.25" (26 cm)  
Blade Width at Widest: 32 mm



## Superior Capsular Retractor

Designed to be placed over the superior capsule.

**PRODUCT NO:**

**7140**  
Overall Length: 9.375" (23.8 cm)  
Blade Width at Widest: 19 mm

