

Mueller Style Hip Instruments



Cross bar can be used to help control rotational stability, and can assist with removal.

PRODUCT NO'S:

6865-01 [Flat Blade Osteotome]
Overall Length: 11.125" (28,3 cm)
Osteotome Width: 20 mm

6865-02 [Femoral Head Dislocation Lever]
Overall Length: 11.375" (23,8 cm)
Scoop Dimensions: 25 mm x 57 mm

6865-03 [Narrow Curved Osteotome]
Overall Length: 12" (30,5 cm)
Osteotome Width: 9 mm

6865-05 [Swan Neck Curved Gouge]
Overall Length: 12" (30,5 cm)
Gouge Width: 23 mm

5350-CB [Cross Bar]



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362