

Proximal Tibia Resection Shark Tooth Clamp

Designed by Shara Diers, PA-C

Designed to help grasp and remove the cut proximal portion of the tibia during total and uni knee arthroplasty

The top looping jaw of the standard size clamp features a sharp forked end for grasping, while the small version features a single pointed end.



PRODUCT NO'S:



3651 [Standard]

Overall Length: 7" (17.8 cm)
Platform Width: .79" (2 cm)
Bottom Jaw Length: 1.5" (3.8 cm)

3651-S [Small]

Overall Length: 7" (17.7 cm)
Platform Width: .25" (0.63 cm)
Bottom Jaw Length: .75" (1.9 cm)

MADE EXCLUSIVELY
FOR INNOMED
GERMANY



ISO 13485:2016

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

© 2025 Innomed, Inc.

103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

INNOMED

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362