

Schanz Pin with Zimmer Hall Quick-connect

Designed by Keith Berend, MD

Used to help remove a femoral head during total hip surgery

Partial threaded pin can be used to help remove a femoral head during total hip surgery. Especially helpful in minimally invasive total hip surgery where access to the femoral head is limited. Connects with a Zimmer Hall Quick-connect.



PRODUCT NO:

3687 [Pin]

Overall Length: 8.625" (21,9 cm)

Shaft Length: 7.375" (18,7 cm)

Thread Length: 2.5" (6,4 cm)

Diameter: 4.5 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

INNOMED