

Universal Tibia/Patella Clamp with Measuring Scale

Designed by S. David Stulberg, MD



Designed to be used to remove a tibia wedge, helps in everting the patella, and the measuring scale can help to assess the thickness of the patella and tibia wedges

PRODUCT NO:

3685

Overall Length: 10" (25,4 cm)
Measuring Scale*: 0 to 26 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



*** Accuracy Statement for devices with a measuring function:**

Devices with a measuring function are considered non-critical, non-calibrated, and are intended for reference-only measurements. While it provides a helpful visual guide for relative measurements, it is not calibrated to specific precision standards. Measurements taken with this device should not be considered accurate enough for critical applications requiring strict dimensional tolerances. General manufacturing tolerances of $\pm 1\%$ at each graduation, or 2% of the full scale apply, but no conformance qualifications are provided or implied.

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362