


Lee MIS Bunion Parallel Guide

Designed by Wonyong Lee, MD


Designed to facilitate MIS bunion correction, focusing on achieving the ideal parallel trajectory of two guide pins for 1st metatarsal screw fixation, with a beveled guide hole to help prevent pin slippage off the bone



PRODUCT NO:

3023

Overall Length: 6.375 (17.1 cm)
Short Tube Guide Length: 1.375" (3.5 cm)
Short Tube Hole Size: For 1.6 mm Pins (.062")
Long Tube Guide Length: 2.375" (6 cm)
Long Tube Hole Size: For 1.4 mm Pins (.054")
Distance Between Holes: 6 mm



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362