



Lombardi Leg Positioner

Designed by Adolph V. Lombardi Jr., MD

Designed to hold the leg during total knee surgery, the unrestricted design helps allow for manipulation of the leg

- ▶ The footpiece consists of an open topped boot, separated with four spacers and attached to a boot support plate.
- ▶ The support plate is bent beyond the rear of the boot to help provide support in different positions of flexion/extension.
- ▶ The spacers allow for room to clean between the boot and support plate.



Two (2) Sterile Pads/Wraps are included with each purchase.

Leg Positioner Sterile Protective Pad & Wrap

Disposable, latex-free sterile foam pad and cohesive wrap helps protect patient from pressure sores, abrasions and possible neurological impairment while securing foot into the boot.

PRODUCT NO'S:	
2629-00	[Case of 10 Sets]
2629-L	[1 Set]



Each set includes
1 Pad and 1 Wrap



PRODUCT NO'S:	
2622	[Lombardi Leg Positioner]
Base Length: 13" (33 cm)	
Base Width: 7" tapers to 4" (17,8 x 10,2 cm)	
Overall Height: 12.75" (32,4 cm)	
Replacement Parts:	
2629-00	[Sterile Pad & Wrap – Case of 10 Sets]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362