

# Posterior Glenoid Neck Retractor

*Used during osteotomy of the humeral head and approaches to the glenoid*

- ▶ Designed to allow one finger retraction
- ▶ Contours to allow teeth to fit behind the glenoid, retracting tissue for easy access to the glenoid



**PRODUCT NO:**

T1002

Width: 30 mm

Overall Length: 10" (25,4 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

# INNOMED

**1.800.548.2362**